





CS500 AUTOMATED GLASS COVERSLIPPER

Introducing the Automatic Glass Coverslipper

This easy to use state of the art Automatic Glass Coverslipper standard with the following features:

- 500 slides/hour out put
- Touch Screen display
- Sensor controlled mountant dispensing
- Self-cleaning pipelines to ensure no blockages
- Auto – detector with alarm for broken coverslips
- Automatic and manual loading mode

 364 Mars-Valencia Rd
Mars, PA 16046

 724.369.1604

 sales@medequipsource.com

 www.medequipsource.com

**CONTACT MEDICAL EQUIPMENT SOURCE
FOR A FREE QUOTE**



Fully Automatic Glass Cover Slipper Specifications

Rated Voltage	a.c. 90V -240V \pm 10%
Rated Power	50Hz
Dimensions (L x W x H)	661 x 726x 627 mm, \pm 5%
Slide Input	90 Slides
Slide Output	90 Slides
Cover Slip Capacity	200 Slides
Mounting Media Volume	250 mL