



## DP360 AUTOMATED INTELLIGENT STAINER

---

### Introducing the Automated Slide Stainer

This easy to use state of the art Automated Stainer comes with the following features:

- High speed throughput
- Minimal cross contamination with the use of a hydrophobic slide carrier
- Multiple Programs can operate simultaneously
- High priority feature to place an emergency batch ahead of the queue
- 4 heated reagent positions
- Exhaust fume monitoring for laboratory safety
- High capacity absorption charcoal filter
- Hybrid high speed stepper motor with optimal sensors for precision movements
- Realtime remote monitoring ability via network
- Intelligent program compensation for reagent consumption
- Automatic report generation for quality control and compliance
- Upgrade without modification to accept an automated coverslip machine



364 Mars-Valencia Rd  
Mars, PA 16046



724.369.1604



[sales@medequipsource.com](mailto:sales@medequipsource.com)



[www.medequipsource.com](http://www.medequipsource.com)

**CONTACT MEDICAL EQUIPMENT SOURCE  
FOR A FREE QUOTE**



## Automated Intelligent Slider Stainer

Rated Voltage	a.c. 110/220V ± 10%
Rated Power	50/60Hz
Max Power Consumption	800 VA
Working Temperature	+15 °C to 40 °C
Dimensions (L x W x H)	1066 x 700 x630mm, ±5%
Total Station	26
Slide Capacity	30 – slides/basket
Reagent Container	480 ± 20 mL
Protocol Storage	>1010