

# MES 1000™


---


## Introducing MEDICAL EQUIPMENT SOURCE'S MES-1000™ Cryostat

This easy to use state of the art Cryostat comes standard with the following features:

- Touch Screen
- UV disinfection
- Hibernation Mode – Sleep Mode
- Automatic or Manual Defrost Cycles
- Expedited Drying Mode
- LED Chamber Light
- Stationary or detachable heat extractor
- Quick Freeze Peltier Shelf -40° C
- Up to 32 Specimen Holders
- Accommodates High- or Low-Profile Blades on same knife block
- Anti-Roll Plate
- Tri-Axial Specimen Clamp



 364 Mars-Valencia Rd  
Mars, PA 16046

 724.369.1604

 sales@medequipsource.com

 [www.medequipsource.com](http://www.medequipsource.com)

**CONTACT MEDICAL EQUIPMENT SOURCE  
FOR A FREE QUOTE**

