



# HP300 TISSUE PROCESSOR

## Introducing the HP300 TISSUE PROCESSOR

This easy to use state of the art Tissue Processor comes standard with the following features:

- 300 Cassette Capacity (Optional 150 Capacity)
- LCD Screen for easy programming and system monitoring
- Automatic self-check of all go/stop functions when unit is on
- Specimen protection function – Automatic infusion of clinically appropriate reagent if unit stops functioning in order to protect specimens from damage
- Remote monitoring via facility Wi-Fi of your choosing
- 4 types of cleaning modes: Quick, Standard, Warm Water and Fat Washing
- Reagent management system
- Heated retort and lid to remove any residual waste via vacuum into waste tray



364 Mars-Valencia Rd  
Mars, PA 16046



724.369.1604



[sales@medequipsource.com](mailto:sales@medequipsource.com)



[www.medequipsource.com](http://www.medequipsource.com)

**CONTACT MEDICAL EQUIPMENT  
SOURCE FOR A FREE QUOTE!**

## HP300 TISSUE PROCESSOR SPECIFICATION

DIMENSION (W x D H)	613 x 680 x 1320mm.
Weight	418lbs
Rate Power:	AC220V ±10% 50/60HZ
Operating Temperature	10~40 °C
Storage Temperature	-10 ~65 °C
Storage of Relative Humidity	30~95%
Number of Reagent Bottles	10
Clean Bottles:	3
Condensate Bottles	1
Paraffin Recycle Bottles	1
Number of paraffin bath(s):	3
Paraffin Volume(s):	4.2L
Number of Paraffin Baths(L):	5.6L
Capacity:	300 Cassettes Maximum
Basket Capacity:	2
Impregnation Vacuum:	0 to-60kPa, Adjustable
Impregnation Pressure:	0 to 60 kPa, Adjustable